



ANALYTICAL TESTING SERVICES

2011 Analytical Price Data Sheet



MEAT SPECIATION	Beef Tissue	Pork Tissue	Poultry Tissue	Sheep Tissue	Horse Tissue	Deer Tissue	Turkey Tissue		Price	Catalog No.	
Cooked Meat Species ELISA											
Cooked Meat Beef	1%								\$75.00	500610	
Cooked Meat Pork		1%							\$75.00	500620	
Cooked Meat Poultry			1%						\$75.00	500630	
Cooked Meat Sheep				1%					\$75.00	500640	
Cooked Meat Horse					1%				\$75.00	500650	
Cooked Meat Deer						1%			\$75.00	500660	
Cooked Meat Turkey							1%		\$75.00	500637	
MELISA-TEK SPECIES											
MELISA Ruminant	0.05%			0.05%					\$75.00	500311	
MELISA Pork		0.05%							\$75.00	500321	
MELISA High Sensitivity Extraction	Ruminant meat and bone meals can be detected at 0.1% (w/w) in feeds and ruminant muscle tissue can be detected at 0.01% in solution.								\$100.00	500391	
THERMAL EVALUATION											
Thermal Evaluation - Pork	Ensures cooking temperature of pork in accordance with EU guidelines.								\$75.00	500629	
RAW MEAT SPECIES											
Cow									Price	Catalog No.	
Raw Meat Cow	1%								\$75.00	500510	
Raw Meat Pig		1%							\$75.00	500520	
Raw Meat Poultry			1%						\$75.00	500535	
Raw Meat Sheep				1%					\$75.00	500540	
Raw Meat Horse					1%				\$75.00	500550	
ADULTERATION											
	milk powder (ppb)										
Bovine Rennet Whey	1%								\$200.00	5171BRW	
FOOD PROTEINS & ALLERGENS											
	Detects:				LOD: (ppm)				Price	Catalog No.	
Almond ²	Specific Almond Proteins				0.5 - 5				\$75.00	500836	
Buckwheat ²	Buckwheat Flour Proteins				2.5 - 25				\$75.00	500806	
Casein ²	Bovine Casein				1 - 10				\$75.00	500811	
Crustacean ²	Shrimp Tropomyosin				0.05 - 0.50				\$75.00	500891	
Egg (Cooked/Raw) ²	Ovomucoid & Ovalbumin				1 - 10				\$75.00	500866	
Gliadin ⁴	Gliadin				5 - 80				\$75.00	500828	
Gliadin Competitive ⁴	Gliadin				5 - 80				\$75.00	500829	
Gluten ³	Gliadin and Glutenins				5				\$75.00	500821	
Gluten ⁵ (AOAC Version)	Wheat Omega Gliadin				≥ 150				\$75.00	500825	
Hazelnut ²	Specific Hazelnut Proteins				0.5 - 5				\$75.00	500839	
Lupin ²	Lupin Flour Protein				0.5 - 5				\$75.00	500855	
Milk/bLG ²	Beta-Lactoglobulin				0.1 - 1				\$75.00	500881	
Peanut ²	Ara h1 & Ara h2				1 - 15				\$75.00	500831	
Prolamins ⁶	Wheat Gliadin & Related Prolamins from Rye, Barley, Triticale				3				\$75.00	500823	

¹ May vary with sample matrix. ² ELISA Systems ³ ELISA Technologies ⁴ R-Biopharm ⁵ HAVEN ⁶ Bio Control
FOB Gainesville, Florida, USA – Volume discounts available – Prices subject to change without notice.



ANALYTICAL TESTING SERVICES

2011 Analytical Price Data Sheet



FOOD PROTEINS & ALLERGENS cont.	Detects:				LOD: (ppm)		Price	Catalog No.		
	Tissue/Muscle (ppb)	Milk (ppb)	Egg (ppb)	Honey (ppb)	Feed (ppb)	Urine (ppb)				
Sesame ²	2S-Albumin Sesame Seed Protein				0.5 - 5		\$75.00	500871		
Soy Protein ²	Soy Tylosin Inhibitor & Soy Flour Proteins				2.5 - 25		\$75.00	500844		
	Tissue/Muscle (ppb)	Milk (ppb)	Egg (ppb)	Honey (ppb)	Feed (ppb)	Urine (ppb)	Serum/Plasma (ppb)	Others (ppb)	Price	Catalog No.
BETA-AGONISTS										
Beta-Agonist	0.2				10	0.75	0.25	0.8 retina 0.25 fecs	\$200.00	5061BAG
Clenbuterol	0.1				5	0.5	0.2	0.5 retina 0.2 feces	\$200.00	5071BAGC
Ractopamine	0.1				4	0.17		0.12 liver	\$200.00	5061BACT
FUNGICIDES										
MG	0.2								\$200.00	5161MGR
ANTHELMINTICS										
Ivermectin	3	5				1	1		\$200.00	5141IVER
Moxidectin	1.5	1.2						20 fat	\$200.00	5841MOXI
ANABOLIC STEROIDS										
Clostebol									\$200.00	5081CLOS
Diethylstilbestrol (DES)	0.9	0.5	0.5			0.5		0.9 liver	\$200.00	5081DES
Ethinylestradiol	0.02					0.04			\$200.00	5841ESTR
Methyltestosterone	0.1					1.3			\$200.00	5841MTES
MedroxyProgesterone Acetate								0.5 fat	\$200.00	5131MPA
Nortestosterone						0.5			\$200.00	5081NOR
Progesterone		0.3				0.3			\$200.00	5081PROG
Stanozolol						1			\$200.00	5081STAN
Trenbolone						3	0.2		\$200.00	5081TREN
Zeranol						1			\$200.00	5081ZER
CORTICOSTEROIDS										
Corticosteroid	0.2	0.2			0.6	3		1 liver	\$200.00	5081COR
Methylprednisolone						0.4			\$200.00	5841PRED
Triamcinolone						0.1			\$200.00	5081TRIA
TRANQUILIZERS										
Azaperone-Azaperol	10							10 kidney 10 liver	\$200.00	5841AZA
Carazolol	1							2 liver 3 kidney	\$200.00	5841CARA
Promazine	0.2							0.2 liver 0.3 kidney	\$200.00	5841PROM
COCCIDIOSTATS										
Diclazuril	2		2						\$200.00	5841DCLA
Halofuginone	0.2		0.1						\$200.00	5841HACO
Ionophores	4				20				\$200.00	5111IONO
Nicarbazin	3.1		0.5						\$200.00	5841NICA

¹ May vary with sample matrix. ² ELISA Systems ³ ELISA Technologies ⁴ R-Biopharm ⁵ HAVEN ⁶ Bio Control
FOB Gainesville, Florida, USA – Volume discounts available – Prices subject to change without notice.



ANALYTICAL TESTING SERVICES

2011 Analytical Price Data Sheet



	Tissue/ Muscle (ppb)	Milk (ppb)	Egg (ppb)	Honey (ppb)	Feed (ppb)	Urine (ppb)	Serum/ Plasma (ppb)	Others (ppb)	Price	Catalog No.
SHELLFISH TOXINS										
Domoic Acid	15							shellfish	\$200.00	5841DOMO
Okadaic Acid	50							shellfish	\$200.00	5841OKA
Saxitoxin	5							shellfish	\$200.00	5841SAXI
MYCOTIXINS ELISA										
	Cereals (ppb)	Milk (ppb)	Wine (ppb)	Nut (ppb)	Feed (ppb)	Food product (ppb)		Others (ppb)	Price	Catalog No.
Aflatoxin DI				0.75		0.02		1.5 peanut butter 0.03 oil	\$200.00	5121AFB
Total Aflatoxin	<0.5								\$200.00	5121AFT
Aflatoxin MI		<0.006						<0.010 cheese 0.016 butter	\$200.00	5121AFM
Deoxynivalenol	30				50	30		beer	\$200.00	5121DON
Fumonisin Bi	2					2			\$200.00	5121FUM
Ochratoxin A	1		0.25			1			\$200.00	5121OCHF
T-2 Toxin						20?50			\$200.00	5121TOX
Zearalenone	6					2	2 liver		\$200.00	5121ZON
ANTIBIOTICS										
	Tissue/ Muscle (ppb)	Milk (ppb)	Egg (ppb)	Honey (ppb)	Feed (ppb)	Urine (ppb)	Serum/ Plasma (ppb)	Others (ppb)	Price	Catalog No.
Bacitracin		10				23			\$200.00	5151BAC
Chloramphenicol	0.02	0.1	0.02	0.02	0.1	0.01		0.01 liver	\$200.00	5091CAPF
Chloramphenicol Fast	0.02	0.02	0.02	0.02	0.5	0.02	0.2	0.02 liver	\$200.00	5091CAPF
Dimetridazole	0.3	0.3	0.3				0.3		\$200.00	5091DIME
Multi-AGP					0.5/1.0				\$200.00	
Tylosin	2.5	2.5		2.5	2.5	2.5			\$200.00	5151TYL
Virginiamycin		8			40	40			\$200.00	5151VIG
AMINO GLYCOCIDES										
Gentamycin	10	2	1	2.5	10	4	2		\$200.00	5111GEN
Neomycin	10	18	14	20		38	35		\$200.00	5111NEO
Streptomycin	10	4	2	5		4	2	5 royal jelly	\$200.00	5111STREP
SULFONAMIDES										
Sulfonamides	4	2.5	3	2		5			\$200.00	5101SULM
Sulfonamides II	5	6		5		13			\$200.00	5101SULMII
Sulfadiazine	7	7		12	12	12	4		\$200.00	5101SUDA
Sulfamethazine	3	8				3	1		\$200.00	5101SUL
Sulfaquinoxaline	5	5				10		5 liver	\$200.00	5101SQUI
FLUOROQUINOLONES										
Difloxacin	8	5							\$200.00	5101DIFX
Enrofloxacin	4	6	9			7	2.5		\$200.00	5101ERFX
Flumequine	0.1	13	3.5		10	5.5	1		\$200.00	5101FLUM
Fluoroquinolones	0.3	3	6	2	5	7	2.5		\$200.00	5101FLUQGG
Fluoroquinolones II	0.6			0.1	16	1.5	3	4 shrimp 6 porcine muscle	\$200.00	5101FLUQII

¹ May vary with sample matrix. ² ELISA Systems ³ ELISA Technologies ⁴ R-Biopharm ⁵ HAVEN ⁶ Bio Control
FOB Gainesville, Florida, USA – Volume discounts available – Prices subject to change without notice.



ELISA Technologies, Inc.
 Analytical Testing Services & Test Kits
 A2LA Accredited to ISO 17025 in Meat Speciation and Allergen Testing



ANALYTICAL TESTING SERVICES

2011 Analytical Price Data Sheet

	Tissue/ Muscle (ppb)	Milk (ppb)	Egg (ppb)	Honey (ppb)	Feed (ppb)	Urine (ppb)	Serum/ Plasma (ppb)	Others (ppb)		
NITROFURANS										
AMOZ	0.1	0.1	0.1	0.1		0.1				5091AMOZ
AOZ	0.05	0.05	0.05	0.05		0.05				5091AOZ
AHD	0.2	0.2	0.2	0.2		0.2				5091AHD
SEM	0.2	0.1	0.1	0.2		0.3		0.1 shrimp		5091SEM

¹ May vary with sample matrix. ² ELISA Systems ³ ELISA Technologies ⁴ R-Biopharm ⁵ HAVen ⁶ Bio Control
 FOB Gainesville, Florida, USA – Volume discounts available – Prices subject to change without notice.